Ottec

Ortesis Térmicas Transpirables

en

Dear customer, thank you for the trust placed in EMO product. Choosing an orthosis of the line Ottec you have acquired a quality product and high health medical level. The Ottec orthotics are made Termotejido: extremely breathable thermal insulation that prevents and dissipates excessive sweating. It acts as an analgesic and muscle relaxant, shortening the time of pre-heating, favoring metabolic processes and increasing blood circulation and exchange of oxygen from the blood into the injured tissues. Read these instructions carefully. If in doubt contact your doctor or specialist dealer.

Normative

The medical device complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices and the respective national regulations. Class I medical device, non-sterile, without measurement function.

Description and characteristics

Wrist MQ510 envelope: holds and protects the wrist joint. Velcro. MQ501 wrist immobilizing Ottec: immobilizes the wrist joint. It incorporates a shaped aluminum removable splint. MQ505 Ottec stabilizing Muñequera thumb: trapezometacarpal stabilizes the joint in opposition, inhibiting the pain during scheduled rest periods. Conformable incorporates a removable aluminum splint on the trapeze-metacarpal joint. Ottec MQ515 Wrist palmar thumb abduction strap: immobilizes the wrist joint and trapezometacarpal joint in abduction, inhibiting the pain during scheduled rest periods. It incorporates a shaped aluminum removable splint on the wrist joint and one on the trapezometacarpal. Velcro CD510: Anatomical shape. Adjustable by inelastic textile shooter. It has a removable pad visco-elastic material. RD510, RD521, RD561, RD564, RD552, RD570, RD575: They have accommodation for the patella, avoiding pressure on it and facilitating flexion and extension. Visco-elastic impeller material. RD561, RD564: Wide adjustable Velcro bands to promote stabilization. metal side strips. RD552: incorporate polycentric joints cap side extension. RD570, RD575: polycentric incorporate adjustable lateral joint control of degrees of flexion indicated. TB501: Open over the calf for easy placement. Incorporates an elastic strap fastener allows the realization of a bandage "eight". TB510: incorporates two lateral rigid straps that stabilize the joint in pronation and supination prevents or investment. Infrapatellar SP515 Support: incorporates a

viscoelastic material cylinder exerts pressure on the patellar tendon, minimizing mechanical commitments. FJ531: lumbosacral orthosis. Close and independent lumbar adjustment for better protection and download lumbar. Contiene four spiral rods side.

Indications

The indications listed below are the most important. All orthoses are indicated for thermocompression treatment, prevention, protection, stabilization, reduction of edema, inflammatory processes, posttraumatic and postoperative rehabilitation, improvement of proprioception and protection during sports activities. Wristbands: minor joint injuries, tendinitis, rheumatic pain, rizoartrosis, bursitis, sprains. If incorporating a splint, immobilization and stabilization of the joint. Coderas: muscular and ligamentous injuries (bursitis, fascitis, rheumatic pain), epicondylitis and epitrocleitis. Shoulders: subluxation of the shoulder joint, instabilities, painful processes. Knee: knee arthritis, tendinitis, inflammatory, mild ligamentous sprains, patellar centering patellar tendon injury. Kneepads with side rods: instability, arthritic and osteoarthritic pain bruises. Kneepads with joints: chronic instability, control of flexion and extension, post-traumatic and post-surgical rehabilitation. Anklets: tendinitis, postural control, medio-lateral (sprains) instability, avoid relapses during sports. Infrapatellar Support: Osgood Schlatter, tendinitis patellar tendon. FJ531 lumbosacral orthoses: LBP, facet syndrome, muscular hypotonia, preventive activities that require a continuous overload of the lumbosacral area.

Instructions for use and placement

The product must be adapted by the prosthetist or by the physician individually for each patient. The right choice of size is crucial for maximum treatment efficacy and prolong the life of the garment. Wrist brace: Adapt the straps to the morphology of the patient. Place the brace on the area to the correct height for use. Help of Velcro strips to close and adjust the garment. epicondylitis Band: place with relaxed muscles and then adjust. Knee: The patellar hole must be properly focused on the kneecap. If incorporate rods or joints, they must be parallel to both sides of the joint. then adjust the Velcro strips. Infrapatellar cinch: place between the patella and the tibial tuberosity. Anklets: put the ankle with the ankle in neutral position of 90°. Subsequently, if incorporated, pronator cinch to adjust

the bandage on "eight" or stabilizer. Trunk orthosis: Place brace on the area to the correct height for use. Close the garment on the abdomen and achieve the desired compression. Subsequently tighten and fix the sidebands in the front of the orthosis, regulating their intraabdominal pressure.

Precautions

Follow the manufacturer's directions for proper orthosis placement. Follow the instructions or recommendations of the prescriber or healthcare personnel. For a personalized adaptation, consult a specialist, if there is ongoing pain, it is advisable to contact a healthcare professional. If an allergic reaction should occur remove the product and consult with the prescribing physician. Do not expose the garment to a flame or heat source at very high temperatures. Store at room temperature, preferably in the original packaging. Guarantee within the framework of laws.

Washing and maintenance

Hand wash in warm water (maximum 30 °) and mild soap. During washing, if the garment has Velcro, they must be glued together or keep them closed. Dry at room temperature. Do not store or iron the garment. Contact with ointments, creams, alcohols and solvents may damage the garment.

Quality certificates

All materials used have been tested, comply with applicable European standards are of high quality and offer maximum safety and comfort during use. Óko-Tex Standard 100: Textile Eco-label awarded after the completion of a catalog of criteria for analysis that guarantees the absence of substances harmful to health. Certifies the absence of harmful substances in all materials with which is made of a garment.

Notice

To dispose of the product, use the solid waste container. For the disposal of the container, depending on its material, use the paper / cardboard container or the plastic container. In any case, you must strictly comply with the regulations of the country. Any serious incident related to the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and / or the patient are established.

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