Ottec

Ortesis Térmicas Transpirables

Dear customer, thank you for placing your trust in an EMO product. By choosing an orthosis from the Ottec line, you have purchased a high-quality product with a high medical and healthcare standard. Ottec orthoses are manufactured with thermal fabric: an extremely breathable thermal insulating material that prevents and dissipates excess perspiration. It acts as an analgesic and muscle relaxant, shortening the pre-warming time, promoting metabolic processes and increasing blood circulation and oxygen exchange from the blood to the injured tissues. Please read these instructions carefully. If you have any questions, contact your doctor or your specialised establishment.

Regulations

This medical device complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices, as well as the applicable national regulations.


Description and characteristics

MQ510 Wrap-around wrist support: supports and protects the wrist joint. Velcro fastening. MQ501 Ottec wrist immobiliser: immobilises the wrist joint. Incorporates a shaped, removable aluminium splint. MQ505 Ottec thumb stabilising wrist support: stabilises the trapeziometacarpal joint in opposition, helping to inhibit pain during scheduled rest periods. Incorporates a mouldable, removable aluminium splint over the trapeziometacarpal joint. MQ515 Ortec thumb abduction wrist support with palmar stay: immobilises the wrist joint and the trapeziometacarpal joint in abduction, helping to inhibit pain during scheduled rest periods. Incorporates a shaped, removable aluminium splint over the wrist joint and another over the trapeziometacarpal joint. Velcro fastening. CD510: Anatomical shape. Adjustable by means of an inelastic textile pull strap. Includes a removable viscoelastic pad. RD558 Open stabilising knee support. RD510, RD521, RD561, RD564, RD552, RD570, RD575: Include a patella opening, avoiding pressure on the kneecap and facilitating flexion-extension. Viscoelastic ring. RD561, RD564: Wide adjustable Velcro straps to promote stabilisation. Lateral metal stays. RD552: incorporates lateral polycentric joints with extension stop. RD570, RD575: incorporate adjustable lateral polycentric joints to control the indicated degrees of flexion. TB501: Open over the calf to facilitate fitting. Incorporates an elastic support strap that enables a figure-of-eight bandage. TB510: incorporates two rigid lateral straps that stabilise the joint in pronation and prevent supination and inversion movements. SP515 Infrapatellar support: incorporates a cylinder made of viscoelastic material that

applies pressure to the patellar tendon, minimising its mechanical stress. FJ531: lumbosacral orthosis. Independent lumbar fastening and adjustment for improved lumbar protection and relief. Contains four lateral spiral stays.

Indications

The indications listed below are the most important. All orthoses are indicated for thermocompression treatment, prevention, protection, stabilisation, reduction of oedema, inflammatory processes, post-traumatic and post-surgical rehabilitation, improvement of proprioception and protection during sports activities. Wrist supports: minor joint injuries, tendinitis, rheumatic pain, rhizarthrosis, bursitis, sprains. If they incorporate a splint, they provide immobilisation and stabilisation of the joint. Elbow supports: muscle and ligament injuries, bursitis, fasciitis, rheumatic pain, epicondylitis and epitrochleitis. Shoulder supports: subluxation of the scapulohumeral joint, instabilities, painful conditions. Knee supports: knee arthritis, tendinitis, inflammatory processes, minor ligament sprains, patellar centring, patellar tendon injury. Knee supports with lateral stays: instability, contusions, arthritic and arthrosic pain. Knee supports with joints: chronic instability, flexion-extension control, post-traumatic and post-surgical rehabilitation. Ankle supports: tendinitis, postural control, medial-lateral instability, sprains, prevention of recurrences during sports practice. Infrapatellar support: Osgood-Schlatter disease, patellar tendon tendinitis. FJ531 lumbosacral orthosis: non-specific low back pain, facet syndrome, muscle hypotonia, preventive use in activities involving continuous overload of the lumbosacral area.

Quality certificates

All materials used have been tested, comply with the applicable European regulations, are of high quality and offer maximum safety and comfort during use. Öko-Tex Standard 100 certificate: textile ecological label granted after applying a catalogue of analysis criteria that guarantees the absence of substances harmful to health. It certifies the absence of harmful substances in all materials used to manufacture the garment.

Instructions for use and fitting

The product must be adapted by the orthopaedic technician or by the healthcare professional on an individual basis for each patient. Correct size selection is essential to achieve maximum treatment effectiveness and to prolong the useful life of the garment. Wrist orthoses: adapt the stays to the patient's morphology. Place the orthosis over the area at the correct height for use. Use the Velcro straps to dose and adjust the garment. Epicondylitis band: fit with the muscles relaxed and then adjust. Knee supports: the patellar opening must be correctly centred over the kneecap. If they include stays or joints, these must remain parallel on both sides of the joint. Then adjust the Velcro straps. Infrapatellar strap: place between the kneecap and the tibial tuberosity. Ankle supports: fit the ankle support with the ankle in a neutral 90° position. Afterwards, if included, adjust the pronation strap for a figure-of-eight or stabilising bandage. Trunk orthoses: place the orthosis over the area at the


correct height for use. Close the garment over the abdomen and achieve the desired compression. Then tighten the lateral bands and secure them to the front part of the orthosis, adjusting the intra-abdominal pressure.

Precautions

Follow the manufacturer's instructions for correct fitting of the orthosis. Follow the instructions or recommendations of the prescriber or healthcare personnel. For a personalised fitting, consult a specialist. If pain persists, it is advisable to contact a healthcare professional. If any allergic reaction occurs, remove the product and consult the prescribing doctor. Do not expose the garment to flames or heat sources at very high temperatures. Store at room temperature, preferably in its original packaging. Warranty within the framework of legal provisions.

Washing and maintenance

Hand wash with lukewarm water, maximum 30°C, and neutral soap. During washing, if the garment has Velcro fasteners, they must be attached to each other or kept closed. Dry at room temperature. Do not hang or iron the garment. Contact with ointments, creams, alcohols and solvents may damage the garment.

All materials used have been tested, comply with the applicable European regulations, are of high quality and offer maximum safety and comfort during use. Öko-Tex Standard 100 certificate: textile ecological label granted after applying a catalogue of analysis criteria that guarantees the absence of substances harmful to health. It certifies the absence of harmful substances in all materials used to manufacture the garment.

Notice

To dispose of the product, use the solid waste container. To dispose of the packaging, depending on its material, use the paper/cardboard container or the plastic container. In all cases, the regulations of the country must be strictly complied with. Any serious incident related to the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient are established.

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