Urban

X1 – X7

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INTRODUCTION

Thank you for choosing an Evoluxion Health product. Our manual wheelchairs can meet different needs and different mobility problems. Please do not use the wheelchairs until you have read and understood the instructions. If you do not understand the instructions, please contact healthcare professionals, dealers or the manufacturer's technicians. Failure to do so may result in personal injury and/or property damage.

INTENDED USE

The wheelchair is designed for use in the home, office, public places, buildings, sidewalks, pedestrian paths and firm surfaces. It is intended to assist the user's mobility and should not be used as a ladder or as a means of traction, transporting objects or similar.

INDICATIONS

The device is specifically indicated for people who (due to a wide variety of possible health problems) are unable to stand and/or walk and therefore need an assistive device for some or all daily tasks.

CONTRAINDICATIONS

It is recommended that the user be accompanied by a caregiver or assistance person. The device is not recommended for people with paralysis or other diseases that cause severe body control problems. These persons require special wheelchairs designed to provide greater stability and support for the user's body.

Caution, fingers or other body parts of the user and/or assistant may get caught between the parts of the device when folding, unfolding or adjusting it. Be very careful when adjusting the device. Observe the dimensions of the spaces between the parts of the device to avoid trapping fingers or other body parts.

WHEELCHAIR COMPONENTS

MEASUREMENTS

CR001 – CR007		41cm	43cm	46cm	48cm
Overall length (cm)		98-101,5	98-101,5	98-101,5	98-101,5
Overall height (cm)		90-92	90-92	90-92	90-92
Overall width(cm)		60,5	63	65,5	68
Overall width folded (cm)		33-30,5	33-30,5	33-30,5	33-30,5
Seat width (cm)		41-40,5	43	46-45,5	48
Seat depth (cm)		44-41,5	44-41,5	44-41,5	44-41,5
Back height (cm)		42-43,5	42-43,5	42-43,5	42-43,5
Floor to Seat height (cm)		48-45,5	48-45,5	48-45,5	48-45,5
Weight capacity (Kg)		120	120	120	120
Measurement of front wheel (cm)		20	20	20	20
Measurement of back wheel (cm)		60	60	60	60
Gross weight (Kg)		17,2	17,2	17,2	17,2
Packaging (cm)		79*31*96	79*31*96	79*31*96	79*31*96

Frame	Single cross brace,
Front wheel:	8"x2", PU tire
Rear wheel:	24" wheel, PU tire
Armrest:	Armrest foldable backward, PU padding
Footrest	Swing in and away footrest, leg support
Seat and backrest	Flame -retardant nylon with sponge
Front Fork:	Plastic
Footrest platform	Plastic

INSTALLATION INSTRUCTIONS

- 1) The initial installation of the parts should be carried out by the company's authorized technical personnel.
- 2) In general, wheelchairs are installed when they leave the factory, folded and packed, and pedals are fitted separately.

List of frame bodies and spare parts:

Item	Accessory name	Quantity
1	Folding frame body	1
2	Movable feet (including pedals)	1 pair
3	Wrench 1	1
4	Lower footrest bands	1

- 3) Here are two simple steps to prepare wheelchairs for use :

Step 1: Expand

Hold the hand with both hands so that the left and right frames are slightly apart. When the two sides of the seat cushion are lightly pressed by hand to the level of the seat cushion, the wheelchair is deployed (as shown in Figure 1)

Note: When unfolding, do not pull hard on the left and right frame to avoid damaging the parts.

When pressing the seat cushion down, do not hold your finger on the seat tube to avoid damaging your finger.

Figure 1

Paso 2: When installing the movable feet, first align the positioning hole in the movable foot and the side bracket on the frame. Then lift part A, turn clockwise to adjust. (Figure 2 and 3).

Figure 2

Figure 3

Paso 3: Position the footrest at the desired height. Hold the safety pin and vertically adjust the footrest to the desired height. When the footrest is in the desired position, tighten the bolt. (Figure 4)


Figure 4

Paso 4: Fold down footrest.

Figure 5

Paso 5: Quick release system. Press the button located in the center of the rear wheel. While holding the center button, pull the wheel outward. To replace the wheel, insert it into the center of the rear axle until you hear a click.

Figure 6

Paso 6: Anti-rollover accesorie.

Figure 7

Using a wheelchair

1) Manual use without assistance

Sit steadily in the wheelchair, then hold the steering wheel of the left and right rear wheels with your hands. Both hands push forward at the same time to move the wheelchair forward, release the hands and repeat the above steps when the wheelchair's speed decreases for a certain distance and it is necessary to push the wheelchair again.

Figura 8

2) Brake method

When the wheelchair is stopped, the user has to hold the left and right brake levers with both hands, and then push the brake lever to the self-locking state of the brake at the same time, so that the rear wheel is completely locked to rotate to the idle state and finish the braking process

Figure 9

3) Method of use in case of an obstacle

In case of an obstacle, the attendant will hold the handles with both hands, lower the position of the rear lower tube of the wheelchair with the foot and tilt the wheelchair backwards so that the front wheel is raised and above the obstacle

Figure 10

4) Method of using the folding armrest

To raise the armrest, press the armrest holder with one hand, lift it off the seat tube and flip it over; to return to the initial position, lower the armrest and press lightly.

Armrest
fastening

Seat tube front
support

Figure 12

- 5) Foldable wheelchair transport method
First turn up the left and right pedals and gently lift both ends of the cushion with both hands.

Figure 13

SAFETY PRECAUTIONS

- 1) Before driving a wheelchair alone, practice should be carried out in different road conditions, accompanied by medical personnel, in order to accumulate sufficient experience to ensure safety.
- 2) **Requirements for driving on outdoor pavements.**
This product can normally drive on conventional concrete or asphalt pavement. On roads with unfamiliar terrain or softer soil, reduce speed; avoid driving on wet and soft grassland, as the tires may slip; avoid driving on soft gravel and beaches; if you are not sure of the road condition, try to avoid driving on them.
- 3) **Requirements for driving on a slope**
When using a wheelchair on a sloping surface, do not tip over or abruptly change the direction of the wheelchair, and the emergency brake cannot stop suddenly when the slope is downhill, so as to avoid the risk of forward tipping. Do not ride on a potentially dangerous slope (such as ground covered with ice and snow, grass and chamfers covered by wet leaves, etc.).
- 4) The brake is to lock the movement of the wheelchair in a safe position, not to be used to brake the wheelchair in motion, in particular to be used to brake the wheelchair downhill.
- 5) When getting in or out of the wheelchair, place the pedal upward (folded). Please do not stand on the pedal when getting on and off the bus. Otherwise, it may cause the wheelchair to tip forward and cause an accident.
- 6) Please do not lean forward to pick up objects from the floor with your hands.
- 7) When transferring from the wheelchair to other seats, you should lean backward as far as possible. This prevents the wheelchair from tipping over.
- 8) The wheelchair is capable of supporting a load of 120 kg on a flat road surface, and when the user's body weight exceeds 120 kg, be careful to use or contact the manufacturer.
- 9) When the user pushes the wheelchair by himself/herself, please keep the speed, and the speed should be kept at 3 ~ 5 km/h.

10) Environmental operating temperature conditions: -25 °C ~ + 50 °C; relative humidity: ≤80%; atmospheric pressure: 86 kPa ~ 106 kPa.

MAINTENANCE - SAFETY INSPECTION LIST

During use of the wheelchair, it should be periodically inspected and maintained according to the following procedures:

- Check that the wheelchair walks in a straight line (do not push it to one side, pull it).
- Check that the front and rear wheels are in good condition
- Check that the armrest bracket attachment is secure, reliable and comfortable to remove, and that the screws are securely locked.
- Check that the seat cushions, etc., are not damaged or deformed; if necessary, replace them.
- Check that the armrest cushions are comfortable. -Clean all wheelchair surfaces..
- Check for signs of corrosion on metal parts and replace if corroded or damaged. si hay signos de corrosión en las piezas metálicas y sustitúyalas si están corroídas o dañadas.

Note 1: The cushions, hands and handrails can be cleaned with a clean, slightly damp cloth. Do not use organic solvents for cleaning. If used by patients, manual wheelchairs should be cleaned once a week. If manual wheelchairs are used by infected patients, they should be cleaned and sterilized with disinfectants (such as 75% alcohol).

TRANSPORT AND STORAGE

- 1) This wheelchair cannot be used as a substitute for transport. It is not permitted to ride in a wheelchair when in any motorised means of transport.
- 2) When using the stairs to transfer wheelchairs between floors, no one must be in a wheelchair. When someone is in a wheelchair, the lift should be used to transfer from one floor to another. If you need to use stairs to transport wheelchairs, be sure to remove the occupants and transport them separately. Do not use escalators to transport wheelchairs between floors. Failure to do so may result in serious personal injury and/or property damage.
- 3) The weight of the wheelchair is 15.2 kg, and the centre of gravity of the wheelchair is in the centre. Be skilled in handling so as not to injure yourself. Only safe and inseparable parts may be used as handling focus. Do not use detachable parts during handling, otherwise personal injury and/or damage to the wheelchair may occur.
- 4) Storage
Storage conditions: ambient temperature -40 °C ~ + 65 °C; relative humidity: 80%; pressure 86 kPa ~ 106 kPa.
Do not store the wheelchair in any rain or water environment. Otherwise, it may result in premature rusting of the wheelchair or damage to its surface decorative material, resulting in serious personal injury and/ or damage to the property.
Do not store or use a wheelchair in the vicinity of open fire or flammable items, which will cause serious personal injury and/ or damage to the property.

AFTER – SALE SERVICE

1. After-sales service for this product is provided by Evoluxion Health and authorised distributors. When you need quality assurance services, please contact your distributor first. If you do not receive satisfactory service from them, please contact us directly. Please state: name of distributor, date of purchase, content of defect, product sequence number, etc. Replaceable parts can be ordered or purchased directly from distributors / suppliers.
2. The following cases do not fall within the scope of free maintenance:
 - Wear parts such as backrest cushion, tyre, armrest pad and plastic parts do not fall under the scope of warranty.

- The failure caused by improper use, improper maintenance and improper storage is not required in the manual.
 - Failure caused by self-assembly and disassembly.
 - Due to other accidental or man-made damage.
- 3.Replacement cycle and replacement method of vulnerable parts: It is recommended that vulnerable parts such as seat back pad, tyre, hand pad and plastic parts should be replaced every year.

WARRANTY

This product is guaranteed for three years from the date of purchase by the user. The manufacturer will supply replacement parts for any part that has broken due to manufacturing faults. Please contact the distributor from whom you purchased the product in case you need to make a claim. Faults resulting from accidental damage, misuse, alteration or wear and tear are not covered by this guarantee. Please contact the dealer from whom you purchased the product to activate the warranty.

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